WOODPECKER

The choice of millions of dentists worldwide

U600

Ultrasonic Scaler

Pleasant & Efficient




Intelligent power output
Experience painless treatment

U600 Ultrasonic Scaler

A global professional scaler manufacturer focuses on technological innovation for 30 years.

International high-tech enterprise
30 years innovation
Woodpecker scaler owns **101** national patents **9** national invention patents



Effective and thorough endodontic irrigation





After irrigation, there is no smear layer on surface but a few amount of dentin debris in the root canal.

Root canal irrigation is the most widely used clinical support methods. Its purpose is to remove the substance inside root canal, dissolve diseased tissue, remove pathogenic microorganisms, and remove smear layer. Lubricate the root canal internal wall. During root canal preparation process, prevent the infectious substance from being pushed to the deep part or out of the apical foramen, reducing the risk of instrument fracture in root canal.


Compared with the traditional syringe irrigation, ultrasonic irrigation can better remove bacteria, dentin debris, and smear layer in root canal. And it can well clean special structures such as the lateral branch canal and root canal narrow area. Currently, ultrasonic irrigation is widely used in endodontic therapy in order to receive a better root canal preparation effect.




Sine wave drive

Avoid instantaneous stimulation caused by power surge, ensuring comfortable treatment.

U600 Mild power output



Rectangular power out of other brand scaler



High-quality handpiece

Aluminum shell can be autoclaved,
not easy to cause cracking or deformation.
Latest handpiece sealing technique enables long durability.
User-friendly design brings better hand feeling.



Automatic water supply system

Easy to change water, eliminating redundant steps and increasing the working efficiency.
Solutions such as hydrogen peroxide, sodium hypochlorite, and chlorhexidine can be used.
Clear 50ml scale mark makes it convenient to formulate solutions.



High configuration

G4	G12	P1	P3	P8	E60	E62
×2	×2	×2	×2	×2	×2	×2

Standard configuration

G4	G12	P1	P3	P8	E60	E62
×1	×1	×1	×1	×1	×1	×1

Note: Woodpecker reserves the right of modifying the tip configuration without further notice.

Supragingival scaling





Implant maintenance





Subgingival scaling





Root canal treatment

Endodontic irrigation





For removal of calcification and bad fillings in root canal





For removal of calcification at root canal orifice and top 1/3 of root canal





Root canal retreatment
Removal of foreign body



Special material (E6, E7)
Thin, flexible, good conductivity to ultrasound



Locate the root canal orifice





Diamond-coated ball head
Locate root canal orifice with security, accuracy and efficiency

Cavity preparation



G30,For preparation between pits and fissures



G31,For the preparation of mesiocclusion cavity



G32,For the preparation of distal cavity



Guilin Woodpecker Medical Instrument Co.,Ltd.
Information Industrial Park, Guilin National High-Tech Zone, Guilin, Guangxi, 541004 P. R. China

Woodpecker Wholesale importer in Spain and Portugal.
RUMAR Cedeira S.L.
C. del Obispo Golfín, 11, Bajo A, 28430 Alpedrete, Madrid
Tlfno.: 91 849 01 04
Tlfno. móvil / Whatsapp: +34 626 327 863

RUMAR guarantees the after-sales service of this product through its distributors / technical service partners.
Always in stock! Delivery in 24 hours